Densitometers

DEN-1 and DFN-1R

Compact and efficient benchtop densitometers for measuring turbidity of cell suspensions in a variety of life-science applications.

The densitometers are designed and factory calibrated to measure turbidity in the range of 0.3 to 15.0 McFarland units (DEN-1) 0.00 to >15.00 McFarland units (DEN-1B) with a small standard deviation. If required, they can deliver a wider measurement range (up to 15.00 McFarland units), but with a greater standard deviation.

- Measurement range: 0.3 to 15.0 McFarland units (DEN-1), 0.00 to >15.00 McFarland units (DEN-1B)
- Measurement time: 1 second
- Precision: ±3%
- Standard deviation at 3.00 McFarland units: ±0.1 McF
- User calibration option





Designed for tubes with an outer diameter of 18mm. 16mm tubes can also be accommodated by using the optional tube adaptor D16 (included).

Factory calibrated - retains calibration without power supply. Can be user calibrated with commercial standards or cell suspensions prepared in the laboratory.

DEN-1 and DEN-1B units are calibrated for operation in range 0.3-15.0 McF and 0->15.00 McF, but it is possible to measure turbidity from 0McF to 15McF (note: the standard deviation of the values increases).



DEN-1B densitometer shown

Bright LED display readings, conveniently shown in McFarland units, are clearly visible

Extremely compact design with small footprint and low profile - fits easily into the smallest workspace

Powered by 3 x AA batteries (DEN-1B only) or via external 12v power supply

Applications:

 Life-science - typical applications include determining concentration of cells (bacterial and yeast cells) in the fermentation process, detecting the susceptibility of micro-organisms to antibiotics, identifying micro-organisms with various test systems, and measuring optical density at fixed wavelength

Densitometers - Technical specification		
• = optional	Densitometers	
	DEN-1	DEN-1B
	h: 75mm d: 115mm w: 165mm weight: 0.7kg	h: 75mm d: 115mm w: 165mm weight: 0.7kg
Light source	LED	
Wavelength Inm	565 ±15	
Range (McFarland units)	0.3 to 15.0	0.00 to 15.00
Precision	±3%	
Measurement time sec	1	
Tube diameter, external diameter mm	18	
with D16 adapter mm	16	
Sample volume ml	>2	
Display/display resolution	LED/0.1 McF	LCD/0.01 McF
Operating temperature range	+4 to 40	
Independent power supply	-	3 x AA batteries
External power supply	Input AC 120-230V, 50/60Hz Output DC 12V	
Input voltage V do	12	
Power consumption W	1 (0.08A)	0.1 (0.007A)

Accessories	
DEN MCF STDS set of liquid McFarland standards in 16mm ø glass tubes. 0.5 / 1.0 / 2.0 / 3.0 / 4.0 can be decanted into an alternative tube (shelf-life no longer valid). Requires D16 adaptor, supplied with DEN-1 / DEN-1B.	•
D16 spare tube adaptor for tubes with 16mm outer diameter. (included as standard)	•



Thank you for reading this data sheet.

For pricing or for further information, please contact us at our UK Office, using the details below.

ﷺ UK Office Keison Products,

P.O. Box 2124, Chelmsford, Essex, CM1 3UP, England.

Tel: +44 (0)330 088 0560

Fax: +44 (0)1245 808399

Email: sales@keison.co.uk

Please note - Product designs and specifications are subject to change without notice. The user is responsible for determining the suitability of this product.