


## A new dimension in flexibility

**DR 2800 Spectrophotometer for water analysis** 



# Incomparable flexibility—systematic water analysis with DR 2800

Wastewater, drinking water and process water—the HACH LANGE system stands for optimal water analysis. A unique selection of perfectly harmonised components guarantees individually accurate solutions. Its rugged design and intuitive handling make the DR 2800 the ideal spectrophotometer for economical, efficient and routine analysis.



## **Competence in water analysis**

#### **Photometer**

DR 2800—the specialist for water analysis in the laboratory and in the field.


- → More than 200 preprogrammed tests
- → Automatic test identification
- → 50 user methods
- → Waterproof (IP 41)

#### **Parameters**

All key parameters from Ammonium to Zinc:

- → Sum parameters—COD/BOD/TOC
- → Nutrients
- → Metals
- → Disinfectants
- → Surfactants and many more
- → Including sample preparation
- → Including quality assurance

For wastewater, drinking water and process water



### THE HACH LANGE SYSTEM FOR WATER ANALYSIS

#### **Services**

- → On-site support by the Technical Customer Service
- → Own service technician and maintenance contracts with warranty extension
- → Legal compliance and environmental protection through collection of used reagents
- → Comprehensive web site with downloads and shop
- → Regular customer information by post and email

#### All methods in one instrument

From fast screening tests to standard-comparable analysis:

DRINKING WATER,
WASTEWATER COOLING WATER AND BOILER WATER

CUVETTE TEST ACCUVAC PERMACHEM

Automatic result assurance

Manual result assurance

## Especially robust in the laboratory and in the field

The reliable optical system of the compact DR 2800 ensures precise results. Each wavelength is individually selectable, 50 user methods and up to 500 measured values can be accommodated in the data storage. Thanks to the variety of cell sizes, a wide range of measured parameters can be covered, from trace analysis to high concentrations in electroplating baths. Easily understandable menus and the large touchscreen minimise familiarisation and prevent errors from occurring.

# The ultimate in simplicity— Cuvette tests with laser identification

Just insert the cuvette—everything else is automatic. A laser identifies the cuvette from its barcode. The DR 2800 calls up the correct test data. As the cuvette rotates, it is measured 10 times and any readings that are distorted by soiling on the cuvette are automatically eliminated. The measurement result is shown immediately in mg/l. All without even having to push a button!



Uniquely flexible: DR 2800 and the proven reagents from HACH LANGE with a wide variety of methods and parameters

## **Technical data**

| DR 2800 Spectrophotometer | Model no. LPV422  |
|---------------------------|---|
| Applications              | Municipal and industrial wastewater, drinking water, sur-   |
|                           | face water, water treatment, process water, production      |
|                           | control   |
| Readout modes             | Absorbance (± 3 E), Transmittance (%), Concentration        |
| Wavelength range          | 340 – 900 nm (tungsten lamp)                                |
| Reference beam technology | To compensate for ageing of the lamp and power              |
|                           | fluctuations  |
| Integrated Barcode        | For automatic identification of LANGE Cuvette Tests, with   |
| Reading System (IBR)      | 10 measurements during one rotation and outlier elimination |
| Cells                     | Cell holder for LANGE Cuvette Tests with IBR system         |
|                           | Cell holder for 10-mm, 50-mm and 1-inch rectangular cells   |
|                           | and 1-inch round cells                                      |
| Display                   | High-resolution LCD (320 x 240 pixel),                      |
|                           | switchable backlighting, touchscreen                        |
| Data storage              | Up to 500 measured values                                   |
|                           | 50 user methods   |
|                           | More than 200 preprogrammed HACH and LANGE method           |
| Ports                     | USB port for PC, printer, USB stick                         |
|                           | Optional: Software for direct data transfer to PC           |
| Protection class          | IP 41   |
| Power requirement         | 100 – 120 V; 200 – 240 V; 50/60 Hz; automatic switching;    |
|                           | including plugs for all European power networks             |
| Battery (optional)        | Lithium battery for up to 40 h of operation                 |
| Ambient temperature       | 10 – 40 °C, 0 – 90 % rel. humidity                          |
| Dimensions (WxHxD)        | 22 x 13.7 x 33.2 cm   |
| Weight                    | 4 kg (without hattery)                                      |

Subject to change.

Nitrate TOC
Hardness TNb Ammonium
COD Phosphate Surfactants
Organic acids Metals Chlorine

And many more



Thank you for reading this data sheet.

For pricing or for further information, please contact us at our UK Office, using the details below.

UK Office Keison Products,

P.O. Box 2124, Chelmsford, Essex, CM1 3UP, England.

Tel: +44 (0)330 088 0560

Fax: +44 (0)1245 808399

Email: sales@keison.co.uk

Please note - Product designs and specifications are subject to change without notice. The user is responsible for determining the suitability of this product.