
Heating Mantles
HM50C, HM100C,
HM250C, HM500C,
HM1000C and HM2000C

Instruction Manual
M8175 Version 1.2

HEATING MANTLE

HM50C, HM100C, HM250C, HM500C,
HM1000C and HM2000C

Instructions for use

Introduction

Thank you for purchasing this piece of Stuart equipment. To get the best performance from the equipment, and for your own safety, please read these instructions carefully before use. Before discarding the packaging check that all parts are present and correct.

This equipment is designed to operate under the following conditions:-

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply fluctuation not exceeding 10%
- ❖ Overvoltage category II IEC60364-4-443
- ❖ Pollution degree 2

If the equipment is not used in the manner described in this manual and with accessories other than those recommended by Stuart, the protection provided may be impaired.

Electrical Installation

THIS INSTRUMENT MUST BE EARTHED

Before connection please ensure that the line supply corresponds to that shown on the rating plate.

| Model | Supply requirements | Power |
|---------|---------------------|-------|
| HM50C | 230V, 50/60Hz | 75W |
| HM100C | 230V, 50/60Hz | 100W |
| HM250C | 230V, 50/60Hz | 200W |
| HM500C | 230V, 50/60Hz | 280W |
| HM1000C | 230V, 50/60Hz | 380W |
| HM2000C | 230V, 50/60Hz | 500W |

All models are supplied with two mains leads fitted with IEC plugs for connection to the instrument. One has a U.K. 3 pin plug and the other has a 2 pin "Schuko" plug for connection to the mains. Choose the lead appropriate for your electrical installation and discard the other. Should neither lead be suitable, take the lead with the U.K. plug and replace the plug with a suitable alternative. This involves cutting off the moulded plug, preparing the cable and connecting to the re-wireable plug in accordance with its instructions.

THIS OPERATION SHOULD ONLY BE UNDERTAKEN BY A QUALIFIED ELECTRICIAN

NOTE: Refer to the equipment's rating plate to ensure that the plug and fusing are suitable for the voltage and wattage stated.

The wires in the mains cable are coloured as follows:

BROWN - LIVE
BLUE - NEUTRAL
GREEN/YELLOW - EARTH

The instruments are fitted with an IEC socket at the rear of the instrument for connection of the mains lead. The appropriate mains lead should be connected **BEFORE** connection to the mains supply.

Should the mains lead need replacement, a cable of 1mm² of harmonised code H05W-F connected to an IEC320 plug should be used.

IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN

Safety Advice Before Use

- ❖ Never move or carry the unit while it is in use or connected to the mains supply.
- ❖ Do not touch the heating element or any glass vessel whilst in use.
- ❖ Only use flasks corresponding to the correct capacity of your heating mantle.
- ❖ Do not use the units with flammable liquids or in a hazardous atmosphere.
- ❖ Do not immerse the unit in water or fluids.
- ❖ Always fill the glassware/vessel away from the heating mantle and ensure there is no liquid on the outside of the vessel before placing it in the mantle. If spillage does occur, disconnect from the mains supply and follow the instructions as detailed in Maintenance & Servicing.
- ❖ Do not cover the mantle while in use.
- ❖ Do not block or obstruct the ventilation slots.
- ❖ Do not leave the equipment switched on without a charged flask.
- ❖ Do not position the unit such that it is difficult to disconnect it from the mains by removing the mains cable from the socket at the back.

Hot Surfaces


WARNING: The product may be HOT. Never move or carry the product until it has been disconnected from the electricity supply and allowed to cool.

General Description


The Stuart heating mantles have been specifically designed for heating fluids in round bottomed flasks. Heating control is provided by a built-in solid state simmerstat (50, 100, 250 and 500ml mantles) or an energy regulator (1000 and 2000ml mantle). The maximum element temperature is 450°C. The mantles are housed in an aluminium case with stove paint finish to give good chemical resistance. The units are provided with ventilation slots in the base and around the rim to allow convection cooling. All mantles also have a support rod positioning clamp at the side of the unit.

Controls

Front view

Rear view

Operation

1. Place the unit on a clean, dry, non-combustible, solid work surface with at least 300mm suitable clearance all around from other equipment.
2. With the mains electricity supply switched off, connect the mains lead with IEC plug to the mains input IEC socket. Connect the plug to the mains supply.
3. Place a charged, clean, dry glass vessel of the size indicated on the mantle rating plate label into the unit. Wherever possible, the glass vessel should be supported within the mantle by means of a support rod and clamp.
4. Switch on the mains electricity supply. Adjust the heating control knob to the required setting. Setting 1 is low heat up; setting 10 is maximum heat up.

NOTE: The 'mains power on' indicator will illuminate. The amber 'heating element on' indicator will illuminate/pulsate when the heaters are in operation.

5. When the process is complete, switch the heating control knob to the off position. Disconnect the mains electricity supply.
6. Remove the charged vessel. Handle hot charged vessels with care.

Maintenance & Servicing

WARNING: Ensure the unit is disconnected from the mains electricity supply before attempting maintenance or servicing.

Periodically clean the instrument using a clean cloth. **DO NOT USE SOLVENTS FOR CLEANING ANY PART OF THIS EQUIPMENT.**

Fuse replacement

The mains fuse holder is located at the rear of the unit within the mains input IEC socket. Refer to the table below for the correct fuse type and rating for your model.

| Model | Total heating power (W) | Fuse rating (A) |
|---------|-------------------------|-----------------|
| HM50C | 75 | 0.5 |
| HM100C | 100 | 0.5 |
| HM250C | 200 | 1.25 |
| HM500C | 280 | 2.5 |
| HM1000C | 380 | 2.5 |
| HM2000C | 500 | 2.5 |

Spillage

In the event of spillage or glassware fracture, do not touch the mantle. Disconnect the product from the mains electricity supply and allow the product to cool. Wearing suitable hand protection, carefully remove any pieces of broken glassware. Wipe off all excess liquid from the mantle and surrounding area using an absorbent soft cloth. Drain off any residual fluid retained in the mantle.

In the case of excessive spillage/ flask fracture, invert the mantle and allow it to drain for minimum of one hour. Then proceed with the following drying out procedure. Place the complete mantle, correct way up, in a heated oven at 50°C for a minimum period of 40 hours. Before further use, the mantle must be subjected to electrical safety testing by competent service personnel.

Servicing

Any repairs or replacement of parts MUST be undertaken by suitably qualified personnel.

For a comprehensive list of parts required by service engineers conducting internal repairs, please contact the Service Department, quoting both the model number and serial number:

Email: cpSERVICE@coleparmer.com

Tel: +44 (0)1785 810475

Only spare parts supplied or specified by Cole-Parmer or its agents should be used. Fitting of non-approved parts may affect the performance of the safety features designed into the instrument. If in any doubt, please contact the Service Department of Cole-Parmer or the point of sale.

For any other technical enquiries, please contact the Technical Support Department:

Email: cpTECHSUPPORT@coleparmer.com

Tel: +44 (0)1785 810433.

Spares and Accessories

| Description | Catalogue Number |
|-------------------------------|------------------|
| Mains cable (UK) | AZ9165 |
| Mains cable (EU) | AZ6747 |
| F500mA fuse pack (pack of 10) | AZ9038 |
| F1.25A fuse pack (pack of 10) | AZ9039 |
| F2.5A fuse pack (pack of 10) | AZ9040 |

Warranty

Cole-Parmer warrants this equipment to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of **three (3)** years. In the event of a justified claim, Cole-Parmer will replace any defective component or replace the unit free of charge. This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by installation, adaptation, modification, fitting of non-approved parts or repair by unauthorised personnel.


Cole-Parmer Ltd.
 Beacon Road,
 Stone, Staffordshire,
 ST15 0SA, United Kingdom
 Tel: +44 (0)1785 810475
 Email: cpservice@coleparmer.com
 Web: www.stuart-equipment.com

Technical Specification

| | Specification |
|------------------------------|---|
| Flask capacity: | |
| HM50C | 50ml |
| HM100C | 100ml |
| HM250C | 250ml |
| HM500C | 500ml |
| HM1000C | 1000ml |
| HM2000C | 2000ml |
| Heating element construction | Thermal insulated element wire stitched into a cartridge construction |
| Maximum element temperature | 450°C. Nominal maximum |
| Case construction | Aluminium |
| Thermal insulation | Rockwool mineral insulation |
| Mains input supply voltage | 230V, 50/60Hz |
| IP rating | IPX0 |

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
Cole-Parmer Ltd.
 Beacon Road,
 Stone, Staffordshire,
 ST15 0SA, United Kingdom
 Tel: +44 (0)1785 810475
 Email: cpsservice@coleparmer.com
 Web: www.stuart-equipment.com

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Notes

This product meets the applicable harmonized standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in its vicinity will meet these standards

and so we cannot guarantee that interference will not occur in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, contact the manufacturer.

EU Declaration of Conformity

| | | | |
|---------------------|--|--------------------|------|
| Product | Laboratory Equipment | File Number | P225 |
| Manufacturer | Cole-Parmer Ltd Beacon Road Stone, Staffordshire ST15 0SA United Kingdom | | |

This declaration of conformity is issued under the sole responsibility of the manufacturer

Object of Declaration

Heating Mantles HM Series
(reference the attached list of catalogue numbers)

The object of the declaration described above is in conformity with the relevant Union Harmonisation Legislation:

| | |
|------------------------------|-------------|
| Low Voltage Directive | 2006/95/EC |
| EMC Directive | 2004/108/EC |
| RoHS Directive | 2011/65/EC |

References to the relevant harmonised standards used or references to the other technical specifications in relation to which conformity is declared:

| | |
|-------------------------|---|
| IEC/EN 61010-1:2001 | Safety requirements for electrical equipment for measurement, control and laboratory use. Part 1: General requirements. |
| IEC/EN 61010-2-010:2003 | Particular requirements for laboratory equipment for the heating of materials. |
| IEC/EN 61326-1:2006 | Electrical equipment for measurement, control and laboratory use. EMC requirements. Part 1: General requirements (Class A). |

Signed for and on behalf of the above manufacturer

| | |
|----------------------------------|--------------------------|
| Additional Information | Year of CE Marking: 2012 |
| Place of Issue | Stone, Staffordshire, UK |
| Date of Issue | January 2012. |
| Authorised Representative | Steve Marriott |
| Title | Technical Director |
| Signature | |

Thank you for reading this data sheet.

For pricing or for further information, please contact us at our UK Office, using the details below.

UK Office

Keison Products,

P.O. Box 2124, Chelmsford, Essex, CM1 3UP, England.

Tel: +44 (0)330 088 0560

Fax: +44 (0)1245 808399

Email: sales@keison.co.uk

Please note - Product designs and specifications are subject to change without notice. The user is responsible for determining the suitability of this product.