Operation Manual (EN)

Diaphragm pump 2-headed

Model:

► RE3011C

Contents

1	Important Information.....	4
1.1	General Information	4
1.2	Target Groups	4
1.3	Intended Use.....	4
1.4	Use for an Unauthorized Purpose	4
1.5	Safety Devices	5
1.6	Meaning of the Warning notes	5
1.7	Product Standards, Safety Regulations.....	5
2	Basic Safety Instructions	6
2.1	General Information	6
2.2	Electricity.....	6
2.3	Mechanical Systems	6
2.4	Hazardous Substances	7
2.5	High Temperatures	7
3	Description	8
3.1	Design	8
3.2	Principle of Operation	8
3.3	Gas ballast	8
3.4	Areas of Application	9
3.5	Pump head circuitry	9
3.6	Materials of the medium-affecting pump parts.....	9
3.7	Scope of Delivery	9
4	Technical Data.....	10
4.1	Dimensions	10
4.2	Intake Pressure / Pumping Speed – Diagram	10
4.3	Device Data.....	11
5	Installation and Operation.....	12
5.1	Unpacking	12
5.2	Installation and Connection.....	12
5.3	Operation	12
5.4	Storage.....	12
5.5	Scrap Disposal	12
6	Maintenance and Servicing.....	13
6.1	General Requirements.....	13
6.2	Maintenance Performed by the User	13
6.2.1	Disassembly.....	14
6.2.2	Assembly.....	14
6.2.3	Test	15
6.3	Maintenance by the Manufacturer	15
6.4	Damage Report.....	15
7	Troubleshooting.....	16
8	Spare Parts Overview	17
8.1	Service kit.....	17
8.2	Spare parts view	17
8.2.1	Spare parts list	18

- Damage Report
- Declaration of Conformity

Important Information

1 Important Information

1.1 General Information

The Diaphragm Pump conform to the following directives:

2006 / 42 / EC	Machinery Directive
2014 / 30 / EC	Electromagnetic Compatibility Directive

The CE sign is located on the rating plate. Observe the binding national and local regulations when fitting the pump into installations.

1.2 Target Groups

This Operating Manual is intended for the personnel planning, operating and maintaining Diaphragm Pump.

This group of people includes:

- Designers and fitters of vacuum apparatus,
- Employees working on commercial laboratory and industrial vacuum technology applications and
- Service personnel for diaphragm pumps

The personnel operating and maintaining the diaphragm pump must have the technical competence required to perform the work that has to be done. The user must authorize the operating personnel to do the work that has to be done. The personnel must have read and understood the complete Operating Manual before using the diaphragm pump. The Operating Manual must be kept at the place of use and be available to the personnel when required.

1.3 Intended Use

- The layout of the diaphragm pump must be appropriate for the conditions of use. The user bears the sole responsibility for this.
- The diaphragm pump may only be operated under the conditions stated
 - in the "Technical Data" section,
 - on the type plate, and
 - in the technical specification for the order concerned.
- Diaphragm pumps are approved for extracting, pumping and compressing gases and vapours. If these gases and vapours are toxic or explosive, then the user must observe the currently valid safety regulations for this application. Special types of diaphragm pumps are available for aggressive and explosive gas mixtures.

1.4 Use for an Unauthorized Purpose

It is forbidden to use the pump for applications deviating from the technical data stated on the type plate or the conditions stated in the supply contract, or to operate it with missing or defective protective devices.

1.5 Safety Devices

Measures such as the following are for the safety of the operating personnel:

- electrical connection with a protective conductor (operating mode S1) and an earthing plug,
- Motor protection switch (thermal) and
- "Hot Surface" label on the pump body - warning notice 

The diaphragm pump must not be operated without these elements.

1.6 Meaning of the Warning notes

Take note of the warning notices. They are in the following box:

	CAUTION ! / WARNING !
Hazard which may lead to serious injuries or material damage.	

1.7 Product Standards, Safety Regulations

The Diaphragm Pump meets the following product standards:

DIN EN ISO 12100:2011-03	Safety of machinery - General principles for design - Risk assessment and risk reduction
DIN EN ISO 13857:2008-06	Safety of machinery - Safety distances to prevent hazard zones being reached by upper and lower limbs
DIN EN 1012-2:2011-12	Compressors and vacuum pumps - Safety requirements - Part 2: Vacuum pumps
DIN EN ISO 2151:2009-01	Acoustics - Noise test code for compressors and vacuum pumps - Engineering method (grade 2)
DIN EN 60204-1:2014-10	Safety of machinery - Electrical equipment of machines - Part 1: General requirements
DIN EN 61000-6-2:2011-06 DIN EN 61000-6-4:2011-09	Electromagnetic compatibility (EMC) - Part 6-2: Generic standards - Immunity for industrial environments Part 6-4: Generic standards - Emission standard for industrial environments
DIN EN 61010-1/A1:2015-04	Safety requirements for electrical equipment for measurement, control and laboratory use - Part 1: General requirements
DIN EN 50110-1:2014-02	Operation of electrical installations
Directive 2012/19/EC	Electrical and electronics - old devices (WEEE)
Directive 2011/65/EC	Dangerous materials in electrical and electronics devices (RoHS)
China - RoHS	Environment protection law - China 2007-03

The following additional safety regulations apply in the FR Germany:

BGV A3	Electrical equipment and operating materials
VBG 5	Power-driven machines
BGR 120	Guidelines for laboratories
BGI 798	Hazard assessment in the laboratory
BGG 919 (VBG 16)	Accident prevention regulations for "compressors"
BGR 189 (BGR 195;192;197)	Use of protective working clothes

Observe the standards and regulations applying in your country when you use the diaphragm pump.

Basic Safety Instructions

2 Basic Safety Instructions

2.1 General Information

Warning notices must be observed. Disregarding them may lead to damage to health and property.

The diaphragm pump must be operated by personnel who can detect impending dangers and take action to prevent them from materialising.

The manufacturer or authorized authorised workshops will only service or maintain the diaphragm pump if it is accompanied by a fully completed damage report. Precise information about the contamination (also negative information if necessary) and thorough cleaning of the diaphragm pump are legally binding parts of the contract.

Contaminated diaphragm pumps and their individual parts must be disposed of in accordance with the legal regulations.

The local regulations apply in foreign countries.

2.2 Electricity

The diaphragm pump of operation mode S1 is supplied. When the location of operation mode S1 devices is changed, please note that the testing must be repeated in accordance with DIN EN 0105, DIN EN 0702 and BGV A2.

The local regulations apply in foreign countries.

Please note the following when connecting to the electrical power supply system:

- The electrical power supply system must have a protective connector according to DIN VDE 0100-410 (IEC 60364-4-41).
- The protective connector must not have any breaks.
- The connecting cable must not be damaged.

2.3 Mechanical Systems

Improper use can lead to injuries or material damage. Observe the following instructions:

- Only operate the diaphragm pump with hoses of the specified dimensions.
- The maximum permissible pressure of 1 bar at the suction connection must not be exceeded.
- Hazardous substances must be separated out as far as this is technically possible before they reach the pump.
- External mechanical stresses and vibrations must not be transmitted to the pump. Only use flexible NW 8 laboratory hoses for connecting diaphragm pumps.
- The overpressure generated at the pressure port must not exceed 1 bar.
- The pump must not be used to suck up fluids. Lay the exhaust pipe so that it slopes downwards, so allowing condensate to flow out of the pump. Collect the condensate and dispose of it in an environmentally compatible manner.
- Prevent dyes exuding.
- Maintain a space of least 20 mm between the pump and adjacent parts in order to enable the pump to cool.

CAUTION !

Solid particles in the pumping medium impair the pumping action and can lead to damage. Prevent solid particles penetrating into the pump.

2.4 Hazardous Substances

The operating company bears the responsibility for the use of the diaphragm pump. Hazardous substances in the gases to be pumped can cause personal injuries and property damage. Pay attention to the warning notices for handling hazardous substances. The local regulations apply in foreign countries.

Combustible Gases

Examine before switching on whether that can form gas combustible gas/air mixtures which can be promoted! Consider the regulations of the guideline 1999/92/EC.

Explosive gases

The Diaphragm pump RE3011C is not certified according to ATEX guidelines 94/9/EC.

Aggressive gases

The Diaphragm pump RE3011C is designed for extracting contaminated gases! Especially aggressive gases have to be explicitly checked for material resistance as described in *chapter 3.6* and, if necessary, modified.

Poisonous gases

Use a separator when pumping poisonous or harmful gases. Prevent such substances from leaking out of the appliance or pump. Treat these substances according to the applicable environmental protection regulations.

Test the strength and leak-tightness of the connecting lines and the connected apparatus. Prevent environmental poisons, e.g. mercury, getting into the diaphragm pumps.

Fulfil the requirements, for example:

- German Hazardous Substances Regulation (GefStoffV) of 01. December 2010
- Regulations 2006/121/EC (classification, packaging and identification of hazardous substances),
- Manufacturer's safety data sheets on hazardous substances.

2.5 High Temperatures

The diaphragm pump may heat up as a result of the temperature of the gas being pumped and through compression heat.

Prevent the following maximum permissible temperatures from being exceeded.

- + 40 °C for the environment, and
- + 60 °C for the gas to be pumped.

The motor for single phase alternating current is protected against overload by an integrated motor protection switch.

Description

3 Description

3.1 Design

The diaphragm pump consists of a pump casing and a drive motor. The pump casing contains the drive unit and two pump heads. Both pump head contains a diaphragm and the work valves. Two pair of pump heads are arranged opposite each other. The pump heads are driven via an eccentric shaft with a connecting rod.

Fig. 1 Diaphragm pump RE3011C

3.2 Principle of Operation

Motor, eccentric shaft and connecting rod set the diaphragms in stroke movement. This changes the size of the space between the diaphragms and pump head (pump chamber). Increasing the size of the pump chamber opens the inlet valve while the outlet valve is closed (intake process).

Decreasing the size of the pump chamber ejects the gas through the outlet valve.

The valves are actuated by the gas being pumped. A large proportion of fluid in the diaphragm pump minimizes the pumping efficiency.

3.3 Gas ballast

When condensable vapours are pumped, they may be compressed above the saturated vapour pressure and condense.

Opening the gas ballast valve **(A)** in the suction line of the last pump stage allows air to flow into the pump chamber. This prevents condensation and flushes the pump clear.

Operation leads to increasing the ultimate pressure and the operating temperature.


Fig. 2 Gas ballast valve

3.4 Areas of Application

The Diaphragm Pump is intended for:

- Pumping and compressing neutral and aggressive gases and vapours.
- Generating a vacuum down to an ultimate pressure < 8 mbar.
- Use in physical and chemical laboratories in trade and industry.
- Use for vacuum filtration, vacuum distillation and vacuum drying, and other vacuum technology applications.

3.5 Pump head circuitry

Two-stage:	Both pump heads are connected in series.
-------------------	--

3.6 Materials of the medium-affecting pump parts

Component	Chemical model - RE3011C (resistant to aggressive gases)
Seal	EPDM
Screw fitting / Connecting elements	PVDF, PP
Valve	PEEK
Diaphragm	Elastomer + PTFE Layer
Vacuum hose	PTFE
Connection head / Pump head	PTFE with carbon-fibre reinforcing *)

*) electrically conductive (with manufacturer's certificate of electrical conductivity)
Material resistance to aggressive media see: Publisher Hoppenstedt Publishing (18. September 2007)


3.7 Scope of Delivery

The scope of delivery is specified in the supply contract.

Technical Data

4 Technical Data

4.1 Dimensions

1	Gas ballast valve	
2	Intake port	(Hose nozzle DN 8 for hose inside diameter 8 mm)
3	Pressure port	(Hose nozzle DN 8 for hose inside diameter 8 mm)

Fig. 3 Dimensions

4.2 Intake Pressure / Pumping Speed – Diagram


Fig. 4 Intake Pressure / Pumping Speed – Diagram

4.3 Device Data

Parameter	Unit	Data
Pumping performance 50/60 Hz DIN 28432 at speed of 1500 rpm	m ³ / h	1.0 / 1.1
	l / min	16.7 / 18
Ultimate pressure at speed of 1500 rpm	mbar	< 8
Ultimate pressure with gas ballast at speed of 1500 rpm		18
Max. inlet pressure	bar	1
Max. outlet pressure		1
Pump stages	-	2
Intake and pressure ports	-	Hose nozzle DN 8 for hose inside diameter 8 mm
Ambient temperature	°C	+ 10 to + 40
Max. Operating gas temperature		+ 60
Bearing	-	maintenance-free
Reference surface sound pressure level DIN EN ISO 2151	dB (A)	< 44
Voltage / Frequency (with motor protection switch, switch and cable)	V, Hz	230, 50/60
Power	kW	0.06
Operating mode	-	S 1
Type of protection DIN EN 60529		IP 54
Class of insulation DIN EN 600034-1		F (160°C)
Weight	kg	6.5
Dimensions (W/D/H)	mm	195 / 235 / 145
Order no.	-	RE3011C

The information presented in this material is based on technical data and test results of nominal units. It is believed to be accurate and reliable and is offered as an aid to help in the selection of products.

It is the responsibility of the user to determine the suitability of the product for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

Bibby Scientific Limited does not warrant, guarantee or assume any obligation or liability in connection with this information.

Installation and Operation

5 Installation and Operation

5.1 Unpacking

Carefully unpack the diaphragm pump.

Check the pump for:

- Transport damage,
- Conformity with the specifications of the supply contract (type, electrical supply data),
- Completeness of the delivery.

Please inform us without delay if there are discrepancies between the delivery and the contractually agreed scope of delivery, or if damage is detected.

Please take note of the general terms of business of the manufacturing firm.

In case of a claim under warranty, the device must be returned in packaging that is suitable for protecting it during transport.

5.2 Installation and Connection

1. Set the diaphragm pump on a flat and horizontal surface.
2. Remove the protective caps on the suction and pressure ports.
3. Prepare the connections.
4. Connect the DN 8 vacuum connector to the suction port.
5. Connect the exhaust pipe to the pressure connection.
6. Connect the diaphragm pump to the electrical supply.

5.3 Operation

Observe the basic safety instructions when using the pump.

The diaphragm pump is switched on and off at the operating switch.

5.4 Storage

The pumps are to be stored in a low-dust, interior room within the temperature range from + 5 to + 40 °C and at a relative air humidity < 90%.

Leave the protective elements on the suction and pressure ports. Another equally good protection may be used.

5.5 Scrap Disposal

	CAUTION !
The Diaphragm pump must be disposed of in accordance with the 2012/19/EC guideline and the specific national regulations. Contaminated Diaphragm pumps must be decontaminated according to the laws.	

6 Maintenance and Servicing

6.1 General Requirements

- Check the pump daily for unusual running noises and heat building up on the surface of the pump.
- We recommend changing the diaphragm after 10,000 operating hours. The user may specify that the exchange be made earlier, depending upon the application process.
- Check the electrical and vacuum connections daily.

6.2 Maintenance Performed by the User

WARNING !

Only perform the work that is described here, and that which is permitted to be done by the user.
 All other maintenance and service work may only be performed by the manufacturer or a dealer authorized by him.
 Beware of the pump parts being possibly contaminated by hazardous substances.
 Wear protective clothing if there is contamination.

Scope of permissible work:

- Loosen and remove the hoses
- Open and remove the pump heads
- Inspect the pump chambers, diaphragms and valves
- Deposits in the inside of the pump must be cleaned out
- Change the diaphragms, valves and seals

Tools required:

Tool kit: Order no. RE3011C/1

consists of:

- Pin type face wrench, adjustable, size 3
- Allan key, size 4
- Open spanner, size 17

Maintenance and Servicing

6.2.1 Disassembly

1. Disconnect the power supply and ensure that it cannot be switched on again.
2. Open the screw clamps of the hoses on the pump body with the SW 17 open spanner.
3. Remove four machine screws **(1)** from each connection head with an Allan key, size 4.
4. Lift off the pressure plate **(2)**, the connection head **(3)** and the pump head **(6)**. The valve gasket **(5)** and diaphragm **(8)** are now freely exposed.
5. Loosen the diaphragm **(8)** at the strain washer **(7)** by turning the size 3 pin type face wrench anticlockwise.
6. Clean the valve gasket **(5)** the pump head **(6)** and the diaphragm **(8)** with a soft cloth and acetone.
7. Check that the drive is in good working order.


Fig. 5 Disassembly, Assembly

	WARNING !
Renew defective parts, if necessary! Wear protective gloves! Parts must be renewed at the intervals stated in this Operating Manual or as specified by the user internally! Do not clean with compressed air!	

6.2.2 Assembly (see Fig. 5)

1. Place the pump so that the diaphragm is lying in a horizontal position.
2. Use the size 3 pin-type face wrench to tighten the pressure disc **(9)**, the diaphragm **(8)** and the strain washer **(7)** with the correct torque of 2 - 4 Nm.
3. Bring the connecting rod **(5)** (see Fig. 6) and the diaphragm **(8)** into the central position.
4. Replace the pump head **(6)**.
5. Insert the valve **(5)** and the O-ring **(4)**.
Ensure that they are lying completely flat. Do not insert the burred side facing the sealing surface. Align the connection head flush with the pin.
6. Replace the pressure plate **(2)** and insert and tighten the 4 machine screws with a torque of 3 - 4 Nm.
7. Reattach the hose connections with clamping ring screw fittings.

6.2.3 Test

- Connect a vacuum measuring device to the suction connector and measure the ultimate pressure.
If the device is working properly, then the figure stated in the technical data must be attained within a maximum of one minute.
- The pump must not make any abnormal noises.
- Moving parts must not touch each other.

6.3 Maintenance by the Manufacturer

Repairs and maintenance going beyond the extent of the work described *in chapter 6.2* or reconditioning or modification may only be performed by the manufacturer or authorized workshops.

	WARNING !
The user shall be liable for the consequences of an incorrect damage report or a contaminated pump. The statements in the damage report are legally binding.	

6.4 Damage Report

You find the form of the damage report to the appendix of this manual.

	WARNING !
Incomplete or incorrectly completed damage reports may endanger the service personnel! Give full information in the damage report, in particular regarding a possible contaminating of the sensor.	

Troubleshooting

7 Troubleshooting

Only manufacturing firm and authorized service workshops may work on the diaphragm pump and their accessories during the warranty period.

Trouble	Cause	Remedy	
		by:	with:
Vacuum pump does not start	No power supply	Qualified electrician	Check electrical installation
	Motor defective	Service workshop	Exchange
	Pump body defective		Repair and/or exchange
Vacuum pump does not generate a vacuum or only an inadequate one	Connected apparatus and/or connecting elements leaking	User or Service workshop	Identify and seal the leak, replace the seals and/or hoses if necessary.
	Vacuum pump leaking		Check the hose connections between the pump heads, replace the hoses and/or fittings if necessary.
	Pump head leaking	Service workshop	Repair and/or exchange
	Diaphragm defective	User or Service workshop	Exchange of the diaphragm (see chapter 6.2)
	Valve defective		Exchange of the valve (see chapter 6.2)
	Vacuum pump dirty		General maintenance / cleaning
	Valves dirty		Cleaning condensates and foreign objects out of the valves.
Running noise	Vacuum pump dirty	User or Service workshop	General maintenance / cleaning
Cable(s)	defective and/or brittle	Qualified electrician	Exchange of the cable(s)

8 Spare Parts Overview

The spare parts lists contain all the spare parts and all the information necessary for ordering. When ordering, please quote the description, quantity and serial number!

	CAUTION !
We are not liable for any damage caused by the installation of any parts not supplied by the manufacturer.	

8.1 Service kit

Designation	Order no.
Service kit	RE3011C/2

The Service kit consists of:

Designation	Piece
O-Ring $\varnothing 28 \times 2$	4
O-Ring $\varnothing 12 \times 2$	3
Valve	4
Diaphragm	2

8.2 Spare parts view


Fig. 6 Exploded view

Spare Parts Overview

8.2.1 Spare parts list

Item no.	Designation	Piece
-	Basic pump complete <i>(consisting of position: 1 – 5)</i>	1
1	- Pump casing	1
-	- Drive complete <i>(consisting of position: 2 – 4)</i>	1
2	- Eccentric	1
3	- Piston rod with ball bearing	2
4	- Close tolerance spacer 25x35x1	3
5	- Alternating-current motor	1
	- IEC power cable extension	1
-	Mains connection cable IEC with plug CEE (D)	1
	Mains connection cable IEC with plug BS (UK)	1
6	Pressure washer	2
7	Diaphragm	2
8	Tightening washer	2
9	Pump head 1	2
10	Valve	4
11	O-Ring EPDM Ø 28 x 2	4
12	Connecting head 1	2
13	Pressure plate	2
14	O-Ring EPDM Ø 12 x 2	3
15	Hose nozzle PP, DN 8 - 1/4"	2
16	Gas ballast valve	1
17	Straight threaded joint with seal edge PVDF, 8 - 1/4"	1
18	Vacuum hose PTFE, 8 / 6x1 mm	0.2 m
19	Threaded elbow joint PVDF, 8 - 1/4"	1
20	Handle	1
21	Rubber-metal pad 20 x 24	2
22	Rubber-metal pad 15 x 6	2

	Bibby Scientific Limited Beacon Road, Stone Staffordshire ST15 0SA, UK	Tel: +44 (0)1785 812121 Fax: +44 (0)1785 810405 Email: info@bibby-scientific.com www.bibby-scientific.com	Return Authorization No. RNA:

Company / Facility Street Post Code / Town	<table border="1"> <tr> <td>Order No.:</td> <td>Date of Purchase:</td> </tr> <tr> <td>Type No.:</td> <td>Serial No.:</td> </tr> <tr> <td colspan="2">Type:</td> </tr> </table>	Order No.:	Date of Purchase:	Type No.:	Serial No.:	Type:	
Order No.:	Date of Purchase:						
Type No.:	Serial No.:						
Type:							

1. Description of the fault: (Which fault occurred?) 				
2. Conditions of use: a) Were the pumps / components in use ? ☐ yes ☐ no b) Which oil was used? * _____ c) Field of application (load, process, use): _____				
3. Type of contamination of pumps / components: a) toxic ** ¹⁾ ☐ yes ☐ no b) corrosive ** ¹⁾ ☐ yes ☐ no c) microbiological ** ¹⁾ ☐ yes ☐ no d) explosive ¹⁾ ☐ yes ☐ no e) radioactive ** ¹⁾ ☐ yes ☐ no f) other harmful substances ¹⁾ ☐ yes ☐ no				
* In the case of rotary vane pumps and roots pumps (Important!)		** Pumps / components will only be accepted if proof is provided that they have been properly cleaned. By arrangement, cleaning can be done by Bibby Scientific Limited.		
Type(s) of substances or process-related reaction products which came into contact with the pumps/ components:				
¹⁾ Commercial name / Product name / Manufacturer	¹⁾ Chemical name (possibly formula)	¹⁾ Danger class	¹⁾ Actions to be taken if harmful substances escape	¹⁾ First aid in case of accidents
1.				
2.				
3.				
4.				
Are the pumps / components free of harmful substances? ☐ yes ☐ no				
4. Legally binding declaration: I/ We declare that the statements here are correct and complete. The pumps / components are shipped in accordance with the legal requirements. Hereby we order the repair.				
Contact: _____		Phone: _____		
Department: _____		Fax: _____		
Date: _____				
Legally binding signature: _____		Company Stamp: _____		
The repair can only be made if a completely filled out damage report is received. Otherwise the work may be delayed. Our general terms of delivery and payment apply.				
¹⁾ Answer all questions or mark with ☐ yes or ☐ no. If a box is not filled out, we automatically suppose that there was a contamination and the required measures are charged to you.				

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Membranpumpe 2-köpfig

Modell:

► RE3011C

**France**

Bibby Scientific
Bâtiment le Deltaparc
Silic pars Nord 2
7 rue du Canal, BP 55437 Villepinte
95944 ROISSY Charles de Gaulle,
France
Tel: +33 (0) 1 48 63 78 00
Fax: +33 (0) 1 48 63 78 01
Email: contact@bibby-scientific.com
www.bibby-scientific.fr

North & South America

Bibby Scientific US
t/a Techne Inc.
3 Terri Lane, Suite 10
Burlington, NJ 08016, USA
Toll Free (in NA): 800-225-9243
Tel: +1 609 589 2560
Fax: +1 609 589 2571
Email: labproducts@techneusa.com
www.bibby-scientific.com

UK

Bibby Scientific Limited
Beacon Road, Stone
Staffordshire
ST15 0SA, UK
Tel: +44 (0)1785 812121
Fax: +44 (0)1785 810405
Email: info@bibby-scientific.com
www.bibby-scientific.com

Middle East

Bibby Scientific Middle East Ltd.
PO Box 27842, Engomi 2433
Nicosia
Cyprus
Tel: +357 22 660 423
Fax: +357 22 660 424
Email: sales@bibby-scientific.me
www.bibby-scientific.com

Asia

Bibby Scientific – Singapore
Prudential Tower, Level 26
30 Cecil Street
Singapore 049712
Tel: +65 6631 2976
Fax: +44 (0)1785 810405
Email: info@bibby-scientific.com
www.bibby-scientific.com

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Declaration of Conformity

Vacuum Pump, Model RE3011C

These products comply with the requirements of the EU Directives listed below:

- 2006/42/EG** EC machinery directive
- 2004/108/EG** EC Electromagnetic Compatibility Directive
- 2011/65/EU** Dangerous materials in electrical and electronics devices (RoHS)
- 2012/19/EU** Electrical and electronics - old devices (WEEE)

Compliance with the requirements of these Directives is claimed by meeting the following standards:

DIN EN ISO12100-2010 - Safety of machinery - General principles for design - Risk assessment and risk reduction

DIN EN ISO13857:2008 - Safety of machinery - Safety distances to prevent hazard zones being reached by upper and lower limbs

DIN EN1012-2:1996 - Compressors and vacuum pumps - Safety requirements - part 2: Vacuum pumps

DIN EN ISO2151:2008 - Acoustics - Noise test code for compressors and vacuum pumps – Engineering method (grade 2)

DIN EN60204-1:2006 - Safety of machinery - Electrical equipment of machines - part 1: General requirements

EN61000-6-2:2005 - Electromagnetic compatibility (EMC) - part 6-2: Generic standards - Immunity for industrial environments


EN61000-6-4:2007 - Electromagnetic compatibility (EMC) - part 6-4: Generic standards - Emission standard for industrial environments environments

DIN EN 50110-1 - Operation of electrical installations

DIN EN 61010-1 - Safety requirements for electrical equipment for measurement, control and laboratory use - part 1: General requirements

The above certificates and reports, from independent test houses, are available upon request.

CE mark affixed '15.

Signed:  (Mr S. Marriott)

Date: 6/8/15

Authority: Technical Director
Bibby Scientific Ltd

Thank you for reading this data sheet.

For pricing or for further information, please contact us at our UK Office, using the details below.

UK Office

Keison Products,

P.O. Box 2124, Chelmsford, Essex, CM1 3UP, England.

Tel: +44 (0)330 088 0560

Fax: +44 (0)1245 808399

Email: sales@keison.co.uk

Please note - Product designs and specifications are subject to change without notice. The user is responsible for determining the suitability of this product.