Material test chambers



Material test chambers with individual programming

With a maximum temperature of 300 $^{\circ}$ C and comprehensive programming options, the M series from BINDER is ideally suited for material and aging testing. The high-performance fan ensures fast heating.



Advantages:

- The specialist for demanding heating profiles
- Adjustable high air exchange rate"Made in Germany" quality

Areas of application:



Buildings Materials Industry



Metal Industry / Engineering



Surface Technology



Features	Customer benefits	Characteristics
APT.line™ temperature technology	 Identical test conditions throughout the chamber interior Independent of specimen size and quantity 	 APT.line for maximum precision ™ Uniform circulation even under full load Homogeneous temperature conditions throughout test specimens
Performance characteristics	 Broad range of applications Short warm up times Suitable for numerous industry standards 	 High air exchange rate Adjustable fan speed Large power reserves Programmable ventilator flap Temperature range up to 300 °C
Inner chamber concept	 Maximum occupational safety Easy loading and unloading of specimen material Easy cleaning 	 Inner chamber made of stainless steel Very tight door closure with 2-point door latch Low heat dissipation due to 60 mm insulation Rack with tilt protection No permanent fixtures
Controller	 Convenient HMI (Human Machine Interface) Convenient documentation All measured values read at a glance 	 MCS screen controller with color display and 25 storable programs Integrated digital continuous-line recorder for monitoring limits and alarm function
Quality	 Reliable devices with long service lives Short delivery times Minimal maintenance and operating costs 	 Premium quality made in Germany Highly automated series production High-quality materials, state-of-the-art production technology High standard according to DIN 12880 (27-point measurement)
Accessories and Services	 Flexible solution in terms of size, type and equipment Numerous options for special applications BINDER INDIVIDUAL for customer-specific solutions Worldwide BINDER Service 	Comprehensive product portfolio Size 53 to 720 liters Additional product lines: Drying ovens, safety drying ovens, vacuum drying ovens, climate chambers Voltage variants and certificates (UL) Various options: Door with viewing window, access ports, reinforced design for heavy loads, Data Logger Kits Worldwide service network

Performance characteristics



- Electronically controlled APT.line™ preheating chamber assuring temperature accuracy and reproducible results
- Temperature range from 5 °C above ambient temperature to 300 °C
- MCS controller with 25 storable programs of 100 sections each for a maximum of 500 program segments
- · User-friendly LCD screen
 - · Easy-to-read menu guide
 - Integrated electronic chart recorder
 - · Variety of options for the graphic display of process parameters
 - · Real-time clock
- Adjustable ramp function via program editor
- · Program-controlled ventilation flap
- High air-exchange rate through high-performance fan
- · Adjustable fan speed
- Exhaust duct Ø 50 mm
- Temperature safety device class 2 (DIN 12880) with visual alarm
- Printer and communication interface RS 422 for use with optional GMP/GLP and FDA guideline 21 CFR Part 11 compliant APT-COM™ DataControlSystem software
- Units up to 115 liters are stackable
- · 2 chrome-plated racks included
- · BINDER test confirmation



M 720

•	Exterior dimensions	
	Width (mm)	1235
	Height (incl. casters) (mm)	1695
	Depth (mm)	865
	Plus door handle and connection (mm)	150
	Wall clearance, rear (mm)	100
	Wall clearance, side (mm)	160
	Exhaust duct (outer Ø mm)	52
	Steam space volume (I)	869
	Number of doors (ea.)	2

•	Interior dimensions	
	Width (mm)	1000
	Height (mm)	1200
	Depth (mm)	600
	Interior volume (I)	720
	Racks (number standard/max.)	2 / 15
	Load per rack (kg)	45
	Permitted total load (kg)	120
	Weight (empty) (kg)	203

•	Temperature data	
	Temperature range approx. 5 °C above ambient temperature to (°C)	300
	Temperature variation	
	at 70 °C (± K)	0,7
	at 150 °C (± K)	1,9
	at 300 °C (± K)	4,6
	Temperature fluctuation (± K)	0,3
	Warm-up time 1)	
	to 70 °C (min.)	7
	to 150 °C (min.)	21
	to 250 °C (min.)	51
	Recovery time after doors were open for 30 sec. 1)	
	at 70 °C (min.)	1
	at 150 °C (min.)	3
	at 300 °C (min.)	5



M 720

•	Electrical data	
	IP protection class acc. to EN 60529	IP 20
	Voltage (± 10%) 50 / 60 Hz (V)	400 / 3N~
	Nominal power (kW)	5,0
	at 70 °C (W)	570
	at 150 °C (W)	1320
	at 300 °C (W)	2600

1) To 98% of the set value

All technical data are specified for units with standard equipment at an ambient temperature of 25 $^{\circ}$ C and a line voltage fluctuation of $\pm 10\%$. The temperature data is determined in accordance to factory standard following DIN 12880, respecting the recommended wall clearances of 10% of the height, width and depth of the inner chamber. All data was determined at 100% of fan speed. All figures are typical average values for series devices. We reserve the right to alter technical specifications at any time.





Access port

With silicone plugs for introducing external measuring instruments into the chamber. Access ports with 10, 30, 50 mm diameters.



Reinforced interior boiler

Includes 2 reinforced shelves for high loads (max. total load 250 kg).



Door lock

Prevents unauthorized access to the process sequences in the chamber.



Calibration certificate & validation

BINDER can significantly reduce the workload in qualifying devices. Nobody knows our devices as well and has as much experience in certifications as we do.



M 720

Access ports with silicone plug, 10, 30, 50, 100 mm	0
Securing elements for additional fastening of racks (1 set of 4 pieces)	0
Anti-slip rubber pads for safe stacking (1 set of 4 pieces)	0
Keyboard lock	0
HEPA fresh-air filter, Class H 14 (according to EN 1822, min. 99.995% for 0.1-0.3 μm particles)	0
Measurement of air change rate according to ASTM D5374, with definition and protocol according to ambient temperature (measurement in factory)	0
Analog temperature output 4 - 20 mA with 6-pin DIN socket (output not adjustable)	0
Additional PT 100 temperature sensor, flexibly installed with external connection, includes LEMO connector (3-pin)	0
Additional measuring channel for digital display of specimen temperature with flexible PT 100 temperature sensor, measured data recorded via RS 422 interface	0
Temperature measurement acc. to DIN 12880 (27 measuring points) at 150 $^{\circ}$ C or at specified temperature with measuring protocol and certificate	0
Factory calibration certificate. Measurement in center of chamber at 150 °C (302 °F) or at specified testing temperature	0
Extension to factory calibration certificate. Each additional measurement at additional measuring point or temperature	0
Data Logger Kit T 350: For continuous temperature recording of 0 °C to 350 °C. Kit includes 1 data logger, PT 100 sensor with 2 m Teflon extension cable and 1 fixture for mounting to the BINDER unit	0
Data Logger converter cable RS 232 to USB 2.0	0
Data Logger Software: Configuration and evaluation software for all BINDER Data Logger Kits, incl. data cable	0
Rack, chrome-plated	0
Rack, stainless steel	0
Reinforced rack, stainless steel, with 1 set of fasteners (4 pieces), max. load 70 kg	0
Shelf, perforated, stainless steel	0
Locking door handle with key	0
Door gasket, FKM (Viton)	0
2 doors each with window 470 x 290 mm and interior lighting (30 W)	0
Reinforced inner chamber, including 2 reinforced racks (maximum total load 250 kg, max. load per rack 70 kg)	0



Thank you for reading this data sheet.

For pricing or for further information, please contact us at our UK Office, using the details below.

UK Office Keison Products,

P.O. Box 2124, Chelmsford, Essex, CM1 3UP, England.

Tel: +44 (0)330 088 0560

Fax: +44 (0)1245 808399

Email: sales@keison.co.uk

Please note - Product designs and specifications are subject to change without notice. The user is responsible for determining the suitability of this product.